# *Brugermanual for AZM DNU*

AZM: Forsyningscentral for medicinsk vakuum i.h.t. EN/ISO 7396-1, type 802T14930

# Specifikationer

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### Declaration of Contamination of Compressors, Vacuum Pumps and Components

The repair and / or servicing of compressors, vacuum pumps and components will be carried out only if a correctly completed declaration has been submitted. **Non-completion will result in delay.** The manufacturer can refuse to accept any equipment without a declaration. **A separate declaration has to be completed for every single component.**

This declaration may be completed and signed only by authorised and qualified staff.

Customer/Dep./Institute : _____	Reason for returning <input checked="" type="checkbox"/> applicable please mark		
Address : _____	<input type="checkbox"/> repair	<input type="checkbox"/> chargeable	<input type="checkbox"/> warranty
Person to contact: _____	<input type="checkbox"/> exchange	<input type="checkbox"/> chargeable	<input type="checkbox"/> warranty
Phone : _____ Fax: _____	restoring goods because of following reason:		
Order number of customer: _____	<input type="checkbox"/> rent/loan	<input type="checkbox"/> for credit	<input type="checkbox"/> against exchange
	<input type="checkbox"/> exchange already received/arranged		
<b>A. Description of the equipment (machine or component)</b>		Ancillary equipment	
Type : _____	_____		
Catalogue number: _____	_____		
Serial number: _____	_____		
Type of oil used : _____	_____		
<b>B. Condition of the equipment</b>			
	No	Yes	No
1. Has the equipment been used	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
2. Drained (Product/service fluid)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3. All openings sealed airtight	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
4. Purged	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
If yes which cleaning agent and which method of cleaning	_____		
	<b>Contamination :</b>		
	No	Yes	
	<input type="checkbox"/>	<input type="checkbox"/>	toxic
	<input type="checkbox"/>	<input type="checkbox"/>	corrosive
	<input type="checkbox"/>	<input type="checkbox"/>	microbiological
	<input type="checkbox"/>	<input type="checkbox"/>	explosive
	<input type="checkbox"/>	<input type="checkbox"/>	radioactive
	<input type="checkbox"/>	<input type="checkbox"/>	other harmful substances
<b>C. Description of processed substances (Please fill in absolutely)</b>			
1. What substances have come into contact with the equipment :			
Trade name and / or chemical term of service fluids and substances processed, properties of the substances According to safety data sheet (e.g. toxic, inflammable, corrosive, radioactive)			
X Tradename: _____	Chemical name: _____	Remains :	
a) _____		<input type="checkbox"/> Yes	<input type="checkbox"/> No
b) _____		<input type="checkbox"/> Yes	<input type="checkbox"/> No
c) _____		<input type="checkbox"/> Yes	<input type="checkbox"/> No
d) _____		<input type="checkbox"/> Yes	<input type="checkbox"/> No
1. Are these substances harmful ? <input type="checkbox"/> <input type="checkbox"/>			
2. Dangerous decomposition products when thermally loaded Which : _____			
Components contaminated by microbiological, explosive or radioactive products will not be accepted without written evidence of decontamination.			
<b>D. Legally binding declaration</b>			
I / we hereby declare that the information supplied on this form is accurate and sufficient to judge any contamination level.			
Name of authorised person (block letters) : _____			
_____	_____	_____	
date	signatur of authorised person	firm stamp	

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# Styring/PLC


The screenshot shows the GRANZOW control interface. At the top, the logo 'GRANZOW<sup>AS</sup>' is displayed in green, with 'ENERTRONIC - DELTA-P' and 'DIN SIKKERHED' below it. A small green icon with a white arrow pointing down is visible at the bottom center. A black line points from the text box on the right to the 'DIN SIKKERHED' text in the screenshot.

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

# Styring/PLC


The screenshot shows a control panel interface for 'AUT. SKIFT PUMPER'. It features a top-left corner with the number '7'. Below this, the text 'TID FOR AUT. SKIFT:' is followed by a green rectangular display showing '99:99'. Underneath, 'DAG FOR AUTO. SKIFT:' is followed by a larger green rectangular display showing 'Days of the week(CW)'. At the bottom of the panel are four navigation buttons: a left arrow, an up arrow, a down arrow, and a right arrow.

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# Styring/PLC

The screenshot shows a control panel interface with a grid background. On the left, a grey panel contains the following elements:


- Top left: **13**
- Top center: **JUSTER GRÆNSE FOR TRYKFALD FILTER**
- Below that: **MAX-TRYKFALD.** followed by a green button labeled **9999 MBar.**
- Below that: **FILTER 1:** and **FILTER 2:**
- Below that: **FØR:** followed by two buttons, each labeled **9999 MBar.**
- Below that: **EFTER:** followed by two buttons, each labeled **9999 MBar.**
- Below that: **INGEN TRYKFALD FILTER** (twice)
- Bottom: Four navigation buttons (left, up, down, right) with arrows.

A white callout box with a black border is positioned to the right of the panel. It contains the text: **Her justeres alarm grænsen for trykfald over filter.** A black line points from the top-right corner of this box to the green **9999 MBar.** button in the panel.

# Styring/PLC

# PI diagram AZM

<p><b>GRANZOW<sup>AS</sup></b>                  E.JBY INDUSTRIVEJ 26 2600 GLOSTRUP 43 202600                  www.granzow.dk info@granzow.dk</p>		Målforh.	Init.	BT
		Dato	21.11.13	
Skejby Sygehus AZM medicinsk vakuum PI-Diagram		Rettet	20.01.14	
		Ordrenr.		
		TEGNINGSNR.		
		14930		